

HumanMilk Workflow


Utilizing Tangent's T13 Medical Grade Tablet

REQUEST EVAL

For use with your HumanMilk Software



Overview



Mobility

The mobility of the T13 tablet gives clinicians the ability to move along their workflow with ease.



Infection-Control


The T13 tablet comes with a hospital-grade, antimicrobial enclosure, allowing clinicians to clean the device with medical-grade cleaners.



Cerner Validated

As a Cerner Validated device, the T13 tablet ensures compatibility with Cerner Millennium applications.

HumanMilk Workflow



Take the Next Step

REQUEST EVAL